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

Bent backend design on the opposite side of needle tip aiming for preventing occurrence of "coring" - tearing off rubber seal on insulin

The diagram illustrates the 'Bent backend' design for insulin needles. It shows a cross-section of a needle inserted into a vial stopper. A circular inset provides a magnified view of the needle tip and the rubber seal. The original design is shown with a straight needle, while the bent backend design is shown with a curved needle. The bent backend design is shown to reduce the occurrence of tearing off the rubber seal, which is labeled as 'coring'.

Terumo's original design
"bent structure"

Reduce the occurrence of
tearing off rubber

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
Outcome – Better User Experience for Better Compliance/Adherence, which Leads to Better QoL

NANOPASS 34G vs 33G comparison

Overall preference of NANOPASS 34G vs 33G after 1 month (n=151)

Reasons for "prefer to use NANOPASS 34G", multiple answers, n=336


Chart created by TERUMO, based on the data from:
T. Asakura and H. Seino, Clinical usefulness of insulin pen needle "NANOPASS needle II Nanopass 34" – Comparison evaluation with NANOPASS needle Nanopass 33. Japanese Journal of Medicine and Pharmaceutical Science. 2013; 69(1): 137-146. Figure 6, 7

Patient reported outcome about glycemic control

Perceived frequency of hyperglycemia/hypoglycemia evaluation in Diabetes Treatment Satisfaction Questionnaire (DTSQ)

"How often have you felt that your blood glucose levels have been unacceptably high/low recently?"


Chart created by TERUMO, based on the data from:
O. Saito, et al. The usability of 34 gauge needle for pen injections in patients using insulin/GLP-1 analog. Posters / Diabetes Research and Clinical Practice 106S1 (2014) S224-PO346

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References:

1. Frid A, et al. New injection recommendations for patients with diabetes. Diabetes Metab. 2010 Sep;36 Suppl 2:S3-18.
2. Tubiana-Rufi N, et al. Short needles (8 mm) reduce the risk of intramuscular injections in children with type 1 diabetes. Diabetes Care. 1999 Oct;22(10):1621-5.
3. Hirsch LJ, et al. The Injection Technique Factor: What You Don't Know or Teach Can Make a Difference. Clin Diabetes. 2019 Jul;37(3):227-233.
4. Birkebaek NH, et al. A 4-mm needle reduces the risk of intramuscular injections without increasing backflow to skin surface in lean diabetic children and adults. Diabetes Care. 2008 Sep;31(9):e65.
5. Polak M, et al. Subcutaneous or intramuscular injections of insulin in children. Are we injecting where we think we are? Diabetes Care. 1996 Dec;19(12):1434-6.
6. Hofman PL, et al. An angled insertion technique using 6-mm needles markedly reduces the risk of intramuscular injections in children and adolescents. Diabet Med. 2007 Dec;24(12):1400-5.
7. Derraik JG, et al. Effects of age, gender, BMI, and anatomical site on skin thickness in children and adults with diabetes. PLoS One. 2014 Jan 21;9(1):e86637.
8. Kalra S, et al. Optimising Insulin Injection Techniques to Improve Diabetes Outcomes. Diabetes Ther. 2023 Nov;14(11):1785-1799.
9. Klonoff DC, et al. Advance Insulin Injection Technique and Education With FITTER Forward Expert Recommendations. Mayo Clin Proc. 2025 Apr;100(4):682-699.
10. Saito D, et al. The Usability of 34 gauge needle for pen injectors in patients using insulin/GLP-1 analog. Posters/Diabetes Research and Clinical Practice. 106S1 (2014); S224-PO346.
11. Asakura T, et al. Fundamental study of NANOPASS II (34G 4mm) needle for pen injectors. Japanese Journal of Medicine and Pharmaceutical Science. 2013, 69(1):147-155.
12. Asakura T & Seino H. Clinical usefulness of insulin pen needle NANOPASS needle 34" – Comparison with NANOPASS needle 33G. Japanese Journal of Medicine and Pharmaceutical Science. 2013, 69(1):137-145.